




**control-x**  
MEDICAL

cxmed.com


# COBRA

## MOTORIZED MOBILE X-RAY

Coming in 2025



Fully Integrated Digital  
X-Ray System with  
**Wireless DAEC**



Effortless Travel  
& **Motorized**  
Telescopic **Column**



**Fast Charging**  
& **Multi-Shift**  
Autonomy



Large LCD  
User Interfaces



**Advanced Sensors**  
for Safe Operation



Standard Generator  
& Tube







**In-bay Charging**  
of 2 DR Panels



**Advanced Remote**  
**Diagnostics**

# SPECIFICATION

## INTUITIVE USER INTERFACE



>400kg



5 km/h




LiFePo  
(extended range option)  
or Lead Acid



110 - 230 VAC  
50/60 Hz



32kW or 50kW,  
320 mAs



0,6 / 1,3 mm dual focus,  
rotating anode,  
130 kV, 107 kHU



Maintenance & Error Alerts,  
Intelligent Battery Management

Control-X Medical cPlc. H-2120 Dunakeszi, Röntgen Street 6., Hungary



+36 1 381 0301



sales@cxmed.com



cxmed.com